



BRASSELER
USA®


BRASSELER'S NEXT GENERATION ELECTRIC SYSTEM

FORZA ELXT



FORZA

**SHORT
SERIES**



BY YOUR SIDE

FORZA ELXT

Introducing Brasseler's next-generation electric motor, the Forza ELXT—an electric motor with unsurpassed power, exceptional performance and cutting-edge safety. The Forza ELXT sets the standard for electric motors.




- ▶ Most powerful micromotor at an incredible 4.2 Ncm of torque! Nearly 30% more power than the industry average!
- ▶ Reduced Risk of Patient Burns with the advanced Safety Check Function
- ▶ Low Cost, Hassle-free installation with the unique split-module space saving design
- ▶ Lightweight and compact motor provides optimal balance and comfort

Forza™ ELXT

- LED fiber optics
- Torque: 4.2 Ncm
- Speed range: 100-40,000 rpm
- Motor weight: 59 grams
- 3-year warranty

MOTOR COMPARISON



Set Includes:

- Control Unit
- PC Box
- Micromotor with Cord
- Mounting Bracket for PC Box
- Power Cord
- AC Adapter (120V)




SHORT SERIES

Brasseler's new Short Series attachments are a revolution in the industry, paving the way for the next generation of electric handpieces. With its ultra-miniature head and neck design and state-of-the-art gearing, the Forza Short Series attachments will change the way you think about electrics forever!

- ▶ Smallest head size in its class at only 8.9 mm in diameter!
- ▶ Most lightweight attachment at ONLY 53 grams
- ▶ Shorter length for optimal balance
- ▶ Specialized Diamond-Like-Carbon (DLC) gearing for superior performance and increased durability
- ▶ Industry-leading 3-year warranty

Forza™ M5

- 1:5 speed increasing
- Lightweight, high-strength titanium
- Fiber optic
- Speed range: 10,000-200,000 rpm
- For friction grip burs



Forza™ M1



- 1:1 gear ratio
- Lightweight, high-strength titanium
- Fiber optic
- Speed range: 2,000-40,000 rpm
- For latch type burs

SHORTER LENGTH


SMALLER HEAD & NECK



Forza™ M5



Average Electric Attachment Size



ELECTRIC ATTACHMENT COMPARISON

Product Name	Weight (g)	Length (mm)	Body Material	Head Height (mm)	Head Diameter (mm)	Warranty (Years)
Forza M5	53	84	Titanium	13.4	8.9	3
Competitor A	73	96	Stainless Steel	13.5	10.3	2
Competitor B	90	95	Stainless Steel	15.5	9.8	2

ELECTRIC MOTOR COMPARISON

Product Name	Micromotor Weight (g)	Micromotor + Attachment Weight (g)	Body Material	Torque (Ncm)	Safety Function	Warranty (Years)
Forza ELXT	59	112	Titanium	4.2	Yes	3
Competitor A	90	180	Stainless Steel	3.5	No	3
Competitor B	69	134	Stainless Steel	3.0	No	3

FORZA™ ELXT ACCESSORIES/COMPONENTS

ELXT Micromotor Only*	5027637U0		Mounting Bracket for PC Box*	5026848U0	
ELXT-U Control Unit Only*	5027634U0		Mounting Bracket for Touch Pad (Optional)	5026847U0	
ELXT-CD Motor Cord - 1.8m*	5027636U0		Mounting Plate for Touch Pad (Optional)	5026849U0	
ELXT-CDL Motor Cord - 2.2m (Optional)	5027635U0		AC Adapter (120V)*	5010627U0	
Continental Motor Cord (Optional)	5027677U0		Power Cord*	5010787U0	

*Included in complete set

► Order online at BrasselerUSA.com or call **800.841.4522**.
In Canada call **800.363.3838**.



BRASSELER
USA® BY YOUR SIDE